

TRAUMA ABDOMEN/PELVIS WWO**DX: ABDOMINAL TRAUMA****PREP: NONE**

Localizer(s)	Frontal & Lateral		
Scan Series	Abdomen/Pelvis	Abdomen/Pelvis	Kidneys
Scan Phase	Noncontrast	Portal venous	Delay
Instructions	Inspiration	Inspiration	Inspiration
Contrast Agent	-	300, 350, or 370	-
Volume/Rate	-	80ml (370) - 100ml (300)/3.0-3.5ml per sec	-
Contrast Delay	N/A	80 sec	3-5 min
Start	Above diaphragms	Above diaphragms	Above adrenals
End	Symphysis pubis	Symphysis pubis	Below kidneys
Scan Type	Helical	Helical	Helical
kVp	110-140	110-130	110-130
Target CTDIvol	<10 mGy	<10 mGy	<10 mGy
Max CTDIvol	25 mGy	25 mGy	25 mGy
Max Scan Time	25 sec	25 sec	25 sec
Pitch	0.9-1.5:1	0.9-1.5:1	0.9-1.5:1
DFOV	To body wall	To body wall	To body wall
Send to PACS	AX 1x1 standard	AX 1x1 standard	AX 1x1 standard
	COR 2x2 standard	COR 2x2 standard	COR 2x2 standard
	SAG 2x2 standard	SAG 2x2 standard	SAG 2x2 standard
		SAG 2x2 bone (trauma only)	