

**PANCREATIC MASS WWO*****DX: PAINLESS JAUNDICE, UNEXPLAINED BILIARY DILATION, CHARACTERIZE PANCREATIC MASS SEEN ON OTHER IMAGING******PREP: 500ML WATER PO 15-30 MIN PRIOR + 1 CUP TO FULL STOMACH AT START OF SCAN***

<b>Localizer(s)</b>	<b>Frontal (Lateral if required)</b>		
<b>Scan Series</b>	<b>Abdomen</b>	<b>Abdomen</b>	<b>Abdomen</b>
<b>Scan Phase</b>	Non-contrast	Arterial	Portal venous
<b>Instructions</b>	Inspiration	Inspiration	Inspiration
<b>Contrast Agent</b>	-	350 or 370	-
<b>Volume/Rate</b>	-	100ml/4.0ml per sec	-
<b>Contrast Delay</b>	-	Aortic peak + 10 sec	80 sec
<b>Start</b>	Above diaphragms	Above diaphragms	Above diaphragms
<b>End</b>	Bifurcation	Bifurcation	Bifurcation
<b>Scan Type</b>	Helical	Helical	Helical
<b>kVp</b>	100-110	100-110	100-110
<b>Target CTDIvol</b>	<10.0 mGy	<12.0 mGy	<12.0 mGy
<b>Max CTDIvol</b>	25 mGy	25 mGy	25 mGy
<b>Max Scan Time</b>	25 sec	25 sec	25 sec
<b>Pitch</b>	0.9-1.75:1	0.9-1.5:1	0.9-1.5:1
<b>DFOV</b>	To body wall	To body wall	To body wall
<b>Send to PACS</b>	AX 1x1 standard	AX 1x1 standard	AX 1x1 standard
	COR 2x2 standard	COR 2x2 standard	COR 2x2 standard
	SAG 3x3 standard	SAG 3x3 standard	SAG 3x3 standard