

PROTOCOL: WRIST ARTHROGRAM**COIL:** Dedicated wrist coil (first choice) or hand/wrist coil. If not present, small flex coil wrapped around.**POSITIONING:** Fingers straight. Pad fingers dorsal and volar so they cannot move. Superman or arm at side. If arm at side careful attention to avoid phase wrap. Extra effort to make patient comfortable.**CLINICAL INDICATIONS/HISTORY:** WRIST PAIN, FRACTURE, LIGAMENT TEAR, ARTHRITIS, INTERNAL DERANGEMENT, INJURY, TFCC TEAR**PLOTTING:** Axial FIRST off of coronal and sagittal localizers. True axial to forearm Coronal = sigmoid notch distal radius to dorsal/radial corner of styloid Sagittal = perpendicular to coronal**COVERAGE:** Distal metaphysis radius/ulna to mid-shaft metacarpals**TIPS:** Note that finger motion causes artifactual signal in tendons

SCAN ORDER	PLANE	IMAGE CONTRAST/WEIGHTING	MODE/SEQ TYPE	Averages (NEX)	BW (Hz/pixel)	ETL	TR RANGE	TE RANGE	TI	FLIP ANGLE	SLICE/GAP (mm)	FOV (cm)	Resolution	Phase Axis	Send to PACS	Fat sat
1	AX	T2FS	2D/FSE	1	150-180	8-12	3500-6000	65-75		>130	3.0/0	8-10	256 x 256	RL	FULL SERIES	WEAK
2	AX	PD	2D/FSE	2	150-180	8-12	3500-6000	38-45		>130	3.0/0	8-10	512 X 320	RL	FULL SERIES	NONE
3	COR	PDFS	2D/FSE	1	150-180	8-12	3500-6000	38-45		>130	1.5-2.0	8-10	256 X 256	SI	FULL SERIES	WEAK
4	COR	PD	2D/FSE	2	150-180	8-12	3500-6000	38-45		>130	1.5-2.0	8-10	512 x 320-384	SI	FULL SERIES	NONE
5	COR	T1FS	2D/FSE	1	150-180	2-3	450-650	8-10		>130	1.5-2.0	8-10	256 x 192	SI	FULL SERIES	WEAK
6	SAG	PDFS	2D/FSE	1	150-180	8-12	3500-6000	38-45		>130	3.0/0	8-10	320 X 288	SI	FULL SERIES	WEAK