

SCAN ORDER	PLANE	IMAGE CONTRAS T/ WEIGHTI NG	MODE	PULSE SEQ	TR RANGE	TE RANGE	TI	FLIP ANGLE	THICKNES S/GAP (mm)	FOV (cm)	Max Pixel (mm) Fr x Ph	PHASE AXIS	SEND TO PACS	Max scan time (target)	
					COVERAGE										
1	AX LOC	T1 BILAT CLOSED	2D	FSE	TOP OF SELLA TO MAXILLARY TEETH	350-790	MIN		>130	3.0 / 0.3	16-18	.8X.8	LR	FULL SERIES	3:00
2	SAG	PD CLOSED	2D	FSE	TMJ R THEN L	>2500	28-40		>130	3.0 / 0.3	8-13	.5X.5	AP	FULL SERIES	4:00
3	SAG	T2 CLOSED	2D	FSE	TMJ R THEN L	>3500	70-100		>130	3.0 / 0.3	8-13	.5X.5	AP	FULL SERIES	4:00
4	SAG	PD OPEN	2D	FSE	TMJ R THEN L	>2500	28-40		>130	2.0 / 0	8-13	.5X.5	AP	FULL SERIES	4:00
5	SAG	T2 OPEN	2D	FSE	TMJ R THEN L	>3500	70-100		>130	3.0 / 0.3	8-13	.5X.5	AP	FULL SERIES	4:00
6+C	COR	T1 BILAT CLOSED +C FS	2D	FSE	TMJ FOSSA	350-790	MIN		>130	3.0 / 0.3	8-13	.5X.5	LR	FULL SERIES	3:00
7+C	AX	T1 BILAT CLOSED +C FS	2D	FSE	TOP OF SELLA TO MAXILLARY TEETH ROOTS	350-790	MIN		>130	3.0 / 0.3	16-18	.8X.8	LR	FULL SERIES	3:00